



HOVERJACK®

Lifting • Evacuation • Vehicle Extraction

FACILITATE PATIENT LIFTING AND EVACUATIONS IN THE MOST CHALLENGING ENVIRONMENTS WITH THE HOVERJACK® AIR PATIENT LIFT.

The HoverJack® Air Patient Lift allows caregivers to safely lift and transport patients in a comfortable and dignified manner. In the event of a patient fall, the Hoverjack inflates to gently lift the patient from the floor to bed or stretcher height in a supine position. When considering Fall Prevention or Emergency Preparedness guidelines, the Hoverjack products are critical components that address the safety of the patient and the caregiver.

Promotes Patient Comfort & Dignity

- Supine lifting feature prevents patient from incurring further injuries after a fall
- Air inflated chambers provide a more comfortable lift compared to other methods
- Available in two widths with a weight capacity of 1200 lbs. to accommodate most patients, including Bariatric

Improves Staff Efficiency & Safety

- Air-assisted technology reduces the number of caregivers required to lift a patient when responding to a fall
- Self-sealing inflation valves facilitate ease of use, while the deflation valves allow users to quickly lower the patient if needed

Facilitates Care In & Out of the Facility

- System flexibility allows for use anywhere in the facility including small spaces, where floor lifts or ceiling lifts cannot be utilized
- Individual chamber design allows users to easily transfer patients to surfaces of varying heights
- Teflon® coated bottom facilitates easy movement across all surfaces, including parking lots

The HoverMatt® air transfer mattress can be used with the same air supply for safe and easy lateral patient transfers from the HoverJack to a bed or stretcher.



HOVERJACK®

Simplify Stairwell Evacuations with the Evacuation HoverJack® Device

The Evacuation HoverJack® Device allows caregivers to perform safe and efficient stairwell evacuations while keeping the patient in a supine position. Adjustable air chambers inflate to provide a stable and comfortable surface for patients, while easily navigating doorways and narrow spaces. The Teflon® coated bottom allows caregivers to effectively transport the Evacuation HoverJack device across rough surfaces and down stairwells. With all of the same features as the standard HoverJack, the Evacuation HoverJack has the following added benefits to accommodate emergency evacuations:

- Cross-over chest straps, dual lateral body straps and foot-end pouch provide patient security during evacuation
- Tapered head-end improves maneuverability around corners and down all stairwell designs, including curving, L-shaped and steep, straight stairs
- Inflatable head-end wedge facilitates adjustable positioning to improve patient's breathing and comfort
- Weight capacity of 700 lbs. for stairwell evacuations and 1200 lbs. for vertical lifts
- Personal belongings compartment built into device offers convenient storage



The HoverJack may also be used in combination with the HoverMatt to safely and easily extract a patient from a vehicle.



HT-Air® 1200 Air Supply

The HT-Air® 1200 Air Supply offers variable speed inflation with six airflow options. One device can be used to inflate all HoverMatt and HoverJack models, as well as HoverTech's positioning devices. The slower air flow options offer improved patient comfort and safety.

HoverJack® Cart

The HoverJack® Cart is a rolling storage system that allows caregivers to easily respond in the event of a patient fall. With a standard size top basket and larger bottom basket, this cart holds the HoverJack, Air Supply and a HoverMatt transfer mattress.



HoverJack® Battery Cart

The HoverJack® Battery Transport System offers a durable, ergonomic solution for portability and storage of the Air200G/Air400G Air Supply and any HoverJack device, featuring the added convenience of a built-in, rechargeable battery.

With its spacious, lightweight and easy to steer design, the HoverJack Battery Transport System is an all-in-one resource for Rapid Response Teams.

HOVERTECH ADVANTAGES

- Innovative market leader focused on customer needs
- Safety and efficacy testing for CPR applications
- Equipment acquisition programs including the HT-Air Rental Program and HoverTech Rewards Program
- Professionally trained, account-dedicated sales and customer service professionals
- Safe Patient Handling Health Care Ergonomics Lifting Program (HELPSM) provides consultative services to assist facilities to develop or improve a SPHM program

In Emergency Situations, the Safety Slide Provides a Way for Caregivers to Evacuate Patients Quickly and Safely.

The Safety Slide™ Emergency Evacuation Slide Sheet allows attending staff to easily move non-ambulatory patients off beds, through hallways and down stairwells using the patient's bed mattress as an evacuation vehicle. The Safety Slide cocoons the patient and mattress together providing a safe and comfortable evacuation.

- Flat profile allows device to remain under mattress for near instant mobility
- Teflon-infused polyester provides a smooth sliding surface that allows the patient, mattress and Safety Slide to easily glide off the bed frame and be transported to safety
- Non-skid Toughtek® panel on top surface, securing straps, and end hood keep mattress attached to Safety Slide during evacuation
- Adjustable straps securely hold individuals of varying size



SAFETYSLIDE™

PRODUCT SPECIFICATIONS

HoverJack Air Patient Lift

Topside Material: Nylon Oxford
Underside Material: Teflon® infused polyester
Latex free

AVAILABLE SIZES

32" W x 70" L – HJ3202T

39" W x 70" L – HJ3902T

Compatible with all HoverTech Air Supplies

Evacuation HoverJack

Topside Material: Nylon Oxford
Foot Pocket: Urethane coated Packcloth
Bottom Chamber: Reinforced with urethane coated Packcloth
Underside Material: Teflon® infused polyester

Size: 32" W x 73" L – HJ32EV

Compatible with all HoverTech Air Supplies

Safety Slide

Topside Material: Nylon Packcloth
Underside Material: Teflon® infused polyester
Non-skid Panel: Toughtek® PVC-based polymer

Size: 36' W x 83" L – HSS-103

Compatible with HoverTech Air Supplies

Air200G (800W), Air400G (1100W), and HTAIR1200/2300 (Variable Speed)

Carts

HoverJack Cart – HJC-200

HoverJack Battery Cart – HJBSC-200*

**Please note this item is not compatible with the HTAIR1200 air supply*



4482 Innovation Way
Allentown, PA 18109

800.471.2776

Fax 610.694.9601

www.HoverMatt.com

Info@HoverMatt.com

BRO_HJ_0916