

STREAMLINE PATIENT HANDLING TASKS

The breathable HoverSling® Repositioning Sheet functions as both a HoverMatt SPU and sling to streamline patient handling tasks and reduce the amount of time and money spent using multiple products. With the HoverSling Repositioning Sheet, caregivers can use one product for lateral transfers, boosting, positioning, turning and vertical lifting. The patient also benefits from trusted HoverMatt air technology that reduces skin shear during transfers and repositioning.



HOVERSLING®
Repositioning Sheet

The breathable HoverSling Split-Leg offers all of the benefits of the HoverMatt SPU Split-Leg for lateral transfers, boosting, turning and repositioning. When used as a lift sling for seated transfers, the patented design cradles the patient's limbs evenly and securely for greater patient comfort. This unique product can also be used as a limb lifter or stabilizer for wound care and Labor and Delivery applications.



HOVERSLING®
Split-Leg