



ADMISSIONS WERE NEVER SO SIMPLE

With the **ook**«SNOW»*ALL*, admitting a patient becomes easier. Through this single platform, institutions can satisfy the needs of bariatric, medical-surgical and palliative care settings.

Its unequalled versatility serves very simple goals : to simplify caregivers' lives and to offer patients more space and comfort. The **ook**«SNOW»*ALL* achieves this by combining the right capabilities with incomparable ease of use.

The bed also allows institutions to standardize throughout their units, reducing cost of ownership, storage space and equipment handling needs.

CAREGIVER TESTIMONIAL



"We no longer have to accustom ourselves to various specialty beds - it is now the bed that adapts to each of our patients! "



INFINITE™ TECHNOLOGY ENSURES A PERFECT FIT TO YOUR NEEDS


MAXIMIZE PATIENT SAFETY AND SAVE STAFF TIME

The **ook**«SNOW»*ALL* combines width adjustment, gap management and low height in order to maximize patient safety and increase staff efficiency.

Its INFINITE™ Technology allows custom positioning within minimum and maximum width. Gaps between the headboard and siderails always remain compliant with the latest 60601-2-52 IEC standard.

This one-step width expansion facilitates the fit to patient needs and saves staff time.

The width monitoring system can also ensure the bed maintains a dictated width that suits the patient.



Customizable width expansion system



Software-integrated width monitoring



Maintained 60mm gaps



Overload protection width expansion system

A BED THAT DOES IT ALL

MAXIMIZING HOSPITAL BED POTENTIAL

With its unique design approach, the **ook^snow[»]ALL** transcends diverse hospital environments by exceeding expectations in bariatric settings as well as med-surg and palliative care.

Outstanding versatility and tremendous value are made possible with the combination of great hospital bed features : low height, higher weight capacity and expandable width.

The **ook^snow[»]** platform - which already focuses on fall prevention, infection prevention and ease of use - is brought to an entirely new level.



EXCEPTIONAL VERSATILITY AND INCOMPARABLE EFFICIENCY




1060 LB / 481 KG WEIGHT CAPACITY
Superior weight capacity ensures institutions benefit from exceptional versatility



INFINITE™ WIDTH ADJUSTMENT
Customizable width expansion improves patient comfort and staff efficiency



LOW HEIGHT
11" low height reduces the impact of falls in med-surg environments while facilitating mobility in bariatric settings



MAXIMIZED SAFETY
Adjustable width maintains a gap of under 60 mm between the headboard and siderail
In compliance with standard IEC 60601-2-52

BENEFIT FROM ALL THE CAPABILITIES OF THE **ook**SNOW»ALL



11" LOW HEIGHT

Low height maintains underbed access for lifts and over bed tables while reducing the adverse effects of a fall and facilitating patient mobility.



SUPERIOR WEIGHT CAPACITY

1060 lb weight capacity contributes to the equipment's versatility in meeting patient needs.



INFINITE™ WIDTH ADJUSTMENT

Starting at 35", a customizable width expansion can accommodate mattresses up to 48".



SAFE GAP MONITORING

Software-integrated system monitors bed width to ensure safe gaps are maintained.



BUILT-IN EFFICIENT SCALE

A smart scale provides accurate weight in any bed and patient position.



USER-FRIENDLY CLEANING ARCHITECTURE

An open architecture, seamless smooth surfaces and liftable deck covers work together to facilitate cleaning.



BOOSTLESS™ BACKREST SYSTEM

Backrest moves backward to reduce shearing and patient repositioning needs.



FALL RISK MANAGEMENT

Fall prevention is increased with the right combination of low height, bed exit technology and ergonomic siderail design.



Smart Nightlights™ provide patients and caregivers with a clear view of the surroundings



Power outlets with USB ports allow patients to keep their mobile devices within reach and reduce the likelihood of a fall



Liftable deck covers simplify cleaning and maintenance



Width expansion system allows width adjustment in one simple step



Inform Safety Message Center shares information to staff and family

DIGNITY FOR ALL PATIENTS

We believe that maintaining dignity in a hospital setting should start with providing state-of-the-art equipment that can suit the specific needs of multiple care environments.

That is why we designed a bed that provides the space and comfort that patients deserve while facilitating their autonomy.

With its appealing look, the **ook**SNOW»ALL also contributes to the creation of a compassionate care atmosphere, setting the pace for recovery.



Customizable width expansion allows for more space and comfort



Remote patient control facilitates autonomy



Low height encourages patient mobility



Integrated patient helper promotes autonomy

ALLOWING CLOSENESS WITH A LOVED ONE

A BED FOR TWO

Many of us have experienced difficult times in which we need to provide comfort to a loved one in a hospital setting.

Whether it is after a surgery, during hospitalization or in a palliative setting, care environments are rarely conceived for this much-needed proximity between patients and those close to them.

We decided to design a product that would allow patients to find solace and comfort near their loved ones.



BELOVED TESTIMONIAL



“There are no words to describe that little bit of happiness and comfort for the time I got to spend with my husband. [...] This bed allowed me to be close to him without squatting like in the regular hospital beds.”

- Brianna Miles*, NL, Canada

*Name has been changed to maintain confidentiality



Toll free number : 1 844 409-4030

umanomedical.com